## calibration

## December 19, 2022

```
[1]: # Notebook for the custom calibration
[2]: import numpy as np
     import plotly.graph_objects as go
     from scipy.signal import find_peaks
     from scipy.optimize import curve_fit
     from utils.get_emsa_data import get_counts_and_name
[3]: # helper functions important for the calibration
     # converts channel to keV
     def channel_to_kev(
         value=None, arr=None, dispersion=0.0100283, offset=21.078, use_offset=True
     ):
         if not use_offset:
             offset = 0
         if value is not None:
             return (value - offset) * dispersion
         if arr is not None:
             return (np.array(arr) - offset) * dispersion
         else:
             raise ValueError("No value or array provided to chennel_to_keV(...)")
     # this function makes one gaussian
     def gaussian(xaxis, mu, amp, std):
         return amp * np.exp(-((xaxis - mu) ** 2) / (2 * std**2))
     # this function makes n gaussians and a m=deg order polynomial on xaxis
     def n_gaussians_and_m_order_bg(xaxis, deg, *args):
         poly = args[: deg + 1]
         gauss = args[deg + 1 :]
         n = int(len(gauss) / 3)
         y = np.zeros(len(xaxis))
         for i in range(n):
             y \leftarrow gaussian(xaxis, gauss[3 * i], gauss[3 * i + 1], gauss[3 * i + 2])
```

```
y += np.polyval(poly, xaxis)
    return y
# this function removes the background linearly, used to make an initial bq_{\sqcup}
 \hookrightarrow guess
def remove_linear_bg(y, peaks, pixel_removal=40, offset=10):
    bg = y.copy()
    for p in peaks:
        # if the peak is too close to the edge, just say the bg is 0
        if p < pixel_removal:</pre>
            bg[:p] = 0
            continue
        else:
            bg[p - pixel_removal : p + pixel_removal] = np.nan
    # make linear interpolation between nan values
    bg = np.interp(np.arange(len(bg)), np.flatnonzero(~np.isnan(bg)), bg[~np.
 ⇔isnan(bg)])
    bg[: int(offset)] = 0 # first 10 channels or the offset should be 0
    return bg
# This function makes the mu, amp, std and polynomial coefficients
def make_fit(array, x, deg=12, prominence=0.01, pixel_removal=50, offset=20):
    # fits one array to n gaussians and a m order polynomial
    peaks, _ = find_peaks(array, prominence=prominence)
    std = np.ones(len(peaks))
    amp = np.ones(len(peaks))
    # must define the function to cure_fit, since deg is a variable which is \square
 →not allowed in the function,
    # because it is not a parameter of the function that curve fit is trying to
 \hookrightarrow fit
    def fit_func_peaks_and_bg(x, *args):
        poly = args[: deg + 1]
        gauss = args[deg + 1 :]
        n = int(len(gauss) / 3)
        y = np.zeros(len(x))
        for i in range(n):
            y \leftarrow gaussian(x, gauss[3 * i], gauss[3 * i + 1], gauss[3 * i + 2])
        y += np.polyval(poly, x)
        return y
    # removing the background linearly for an initial quess of the polynomial
    bg = remove_linear_bg(array, peaks, pixel_removal=pixel_removal,__

offset=offset)
    # fitting the polynomial
```

```
poly_init = np.polyfit(x, bg, deg)
    # fitting
    init_vals = list(poly_init)
    for i in range(len(peaks)):
        init_vals += [peaks[i], amp[i], std[i]]
    fit_vals, covar = curve_fit(fit_func_peaks_and_bg, x, array, p0=init_vals)
    return fit_vals, covar
def calculate_calibration(peaks_channel, peaks_keV):
    Calibration of the channel width and offset using two peaks.
    The given channels and keV must correspond.
    Parameters
    _____
    peaks_channel : list
        channel value of peaks, fitted
    peaks_keV : list
        theoretical energy of the two peaks
    Returns
    _____
    float, float
        dispersion, offset
    # figure out the distances between the peaks
    channel_distance = peaks_channel[1] - peaks_channel[0]
    kev_distance = peaks_keV[1] - peaks_keV[0]
    # dispersion = (p1_keV - p0_keV) / (p1_channel - p0_channel)
    dispersion = kev_distance / channel_distance # kev_per_channel
    offset = peaks_channel[0] - peaks_keV[0] / dispersion #_
 \hookrightarrow calibrated_zero_channel
    print(
        f"Calibration: {dispersion:.08f} keV/channel, {offset:.03f} channels_
 ⇔zero offset"
    )
    return dispersion, offset
```

```
[4]: # variables used channels = np.arange(2048)
```

```
ga_La = 1.098
as_Ka = 10.5436
peaks_keV = [ga_La, as_Ka]
initial_dispersion = 0.01
initial_offset = 20
deg = 12
prominence = 0.01
file = "GaAs_30kV"
ga30, _ = get_counts_and_name(file)
spectrum = ga30
```

Calibration: 0.01003001 keV/channel, 21.127 channels zero offset



